

TLR4 Polyclonal Antibody RBITC Conjugated

Catalog No: #C90453RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TLR4 Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	751-835/835
Conjugates	RBITC
Target Name	TLR4
Other Names	TLR4
Accession No.	Uniprot:Q9QUK6Gene ID:7099
Uniprot	Q9QUK6
GeneID	7099
Calculated MW	90
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only