

Nephrin Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90454PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nephrin Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	451-550/1241
Conjugates	PE-Cy5
Target Name	Nephrin
Other Names	Nephrin
Accession No.	Uniprot:O60500Gene ID:4868
Uniprot	O60500
GeneID	4868
Calculated MW	138
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only