

## Sialoadhesin Polyclonal Antibody PE Conjugated

Catalog No: #C90456PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sialoadhesin Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms
Immunogen Type	1551-1600/1709
Conjugates	PE
Target Name	Sialoadhesin
Other Names	Sialoadhesin
Accession No.	Uniprot:Q9BZZ2Gene ID:6614
Uniprot	Q9BZZ2
GeneID	6614
Calculated MW	186
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only