

IL-2 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C90467Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-2 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	31-130/153
Conjugates	Cy3
Target Name	IL-2
Other Names	IL-2
Accession No.	Uniprot:P60568Gene ID:3558
Uniprot	P60568
GeneID	3558
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only