

PIK3 gamma Polyclonal Antibody Biotin Conjugated

Catalog No: #C90468Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PIK3 gamma Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	521-620/1102
Conjugates	Bio
Target Name	PIK3 gamma
Other Names	PIK3 gamma
Accession No.	Uniprot:P48736Gene ID:5294
Uniprot	P48736
GeneID	5294
Calculated MW	121
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only