

NF-L Polyclonal Antibody Biotin Conjugated

Catalog No: #C90470Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-L Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-130/543
Conjugates	Bio
Target Name	NF-L
Other Names	NF-L
Accession No.	Uniprot:P07196Gene ID:4747
Uniprot	P07196
GeneID	4747
Calculated MW	68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only