

## HMGA1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90478Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | HMGA1 Polyclonal Antibody Biotin Conjugated                                     |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | WB IHC-P IHC-F IF   |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 51-107/107  |
| Conjugates         | Bio   |
| Target Name        | HMGA1   |
| Other Names        | HMGA1   |
| Accession No.      | Uniprot:P17096Gene ID:3159  |
| Uniprot            | P17096  |
| GeneID             | 3159  |
| Calculated MW      | 12  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only