

TPOR Polyclonal Antibody APC Conjugated

Catalog No: #C90491APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TPOR Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/635
Conjugates	APC
Target Name	TPOR
Other Names	TPOR
Accession No.	Uniprot:P40238Gene ID:4352
Uniprot	P40238
GeneID	4352
Calculated MW	68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only