

Sarcomeric Alpha Actinin Polyclonal Antibody Biotin Conjugated



Catalog No: #C90492Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sarcomeric Alpha Actinin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	351-450/894
Conjugates	Bio
Target Name	Sarcomeric Alpha Actinin
Other Names	Sarcomeric Alpha Actinin
Accession No.	Uniprot:P35609Gene ID:88
Uniprot	P35609
GeneID	88
Calculated MW	98
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only