

## GLUT2 Polyclonal Antibody APC Conjugated

Catalog No: #C90497APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GLUT2 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/524
Conjugates	APC
Target Name	GLUT2
Other Names	GLUT2
Accession No.	Uniprot:P11168Gene ID:6514
Uniprot	P11168
GeneID	6514
Calculated MW	54
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only