

## HBeAg Polyclonal Antibody Biotin Conjugated

Catalog No: #C90498Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HBeAg Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-F
Species Reactivity	Hepatitis B virus
Immunogen Type	51-212/212
Conjugates	Bio
Target Name	HBeAg
Other Names	HBeAg
Accession No.	Uniprot:Q89714Gene ID:944568
Uniprot	Q89714
GeneID	944568
Calculated MW	12
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-F=1:50-200

Note: This product is for in vitro research use only