

TIMP-2 Polyclonal Antibody APC Conjugated

Catalog No: #C90501APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIMP-2 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/220
Conjugates	APC
Target Name	TIMP-2
Other Names	TIMP-2
Accession No.	Uniprot:P16035Gene ID:7077
Uniprot	P16035
GeneID	7077
Calculated MW	24
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only