

Thrombopoietin Polyclonal Antibody FITC Conjugated

Catalog No: #C90505FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thrombopoietin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/353
Conjugates	FITC
Target Name	Thrombopoietin
Other Names	Thrombopoietin
Accession No.	Uniprot:P40225Gene ID:7066
Uniprot	P40225
GeneID	7066
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only