

Collagen I Polyclonal Antibody PE Conjugated

Catalog No: #C90511PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen I Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1051-1150/1464
Conjugates	PE
Target Name	Collagen I
Other Names	Collagen I
Accession No.	Uniprot:P02452Gene ID:1277
Uniprot	P02452
GeneID	1277
Calculated MW	130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only