

Caspase-1 p20 Polyclonal Antibody APC Conjugated

Catalog No: #C90514APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Caspase-1 p20 Polyclonal Antibody APC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 181-280/402 |
| Conjugates | APC |
| Target Name | Caspase-1 p20 |
| Other Names | Caspase-1 p20 |
| Accession No. | Uniprot:P20466Gene ID:12362 |
| Uniprot | P20466 |
| GeneID | 12362 |
| Calculated MW | 20/46 |
| Formulation | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4 |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only