

NSE Polyclonal Antibody FITC Conjugated

Catalog No: #C90515FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	NSE Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF,FC
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/434
Conjugates	FITC
Target Name	NSE
Other Names	NSE
Accession No.	Uniprot:P09104Gene ID:2026
Uniprot	P09104
GeneID	2026
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:100-500 Flow-Cyt=1:100-500/10⁵ cells/Test ICC=1:100-500 IF=1:100-500

Note: This product is for in vitro research use only