

# Complement C5 beta chain Polyclonal Antibody FITC Conjugated



Catalog No: #C90525FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Complement C5 beta chain Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	51-150/1676
Conjugates	FITC
Target Name	Complement C5 beta chain
Other Names	Complement C5 beta chain
Accession No.	Uniprot:P02748Gene ID:727
Uniprot	P02748
GeneID	727
Calculated MW	72
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only