

TTF1 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C90533Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TTF1 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	121-220/371
Conjugates	Cy3
Target Name	TTF1
Other Names	TTF1
Accession No.	Uniprot:P43699Gene ID:7080
Uniprot	P43699
GeneID	7080
Calculated MW	38
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only