

beta Amyloid 1-16 Polyclonal Antibody HRP Conjugated

Catalog No: #C90536HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	beta Amyloid 1-16 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu
Immunogen Type	1-16/42
Conjugates	HRP
Target Name	beta Amyloid 1-16
Other Names	beta Amyloid 1-16
Accession No.	Uniprot:P05067Gene ID:351
Uniprot	P05067
GeneID	351
Calculated MW	4.3/85
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only