

WASF3 Polyclonal Antibody FITC Conjugated

Catalog No: #C90540FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	WASF3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/502
Conjugates	FITC
Target Name	WASF3
Other Names	WASF3
Accession No.	Uniprot:Q9UPY6Gene ID:10810
Uniprot	Q9UPY6
GeneID	10810
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only