

CD45 Polyclonal Antibody FITC Conjugated

Catalog No: #C90548FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD45 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Rat
Immunogen Type	131-230/1273
Conjugates	FITC
Target Name	CD45
Other Names	CD45
Accession No.	Uniprot:P08575Gene ID:24699
Uniprot	P08575
GeneID	24699
Calculated MW	143
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only