

Myeloperoxidase Polyclonal Antibody FITC Conjugated



Catalog No: #C90552FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myeloperoxidase Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/745
Conjugates	FITC
Target Name	Myeloperoxidase
Other Names	Myeloperoxidase
Accession No.	Uniprot:P05164Gene ID:4353
Uniprot	P05164
GeneID	4353
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only