

EGFRvIII Polyclonal Antibody FITC Conjugated

Catalog No: #C90553FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EGFRvIII Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	321-520/1210
Conjugates	FITC
Target Name	EGFRvIII
Other Names	EGFRvIII
Accession No.	Uniprot:P00533Gene ID:1956
Uniprot	P00533
GeneID	1956
Calculated MW	100/130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only