

# HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated



Catalog No: #C90555PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated                  |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | ICC IF  |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 151-250/365   |
| Conjugates         | PE-Cy5  |
| Target Name        | HLA Class 1 ABC /HLA ABC  |
| Other Names        | HLA Class 1 ABC /HLA ABC  |
| Accession No.      | Uniprot:P30443Gene ID:3105  |
| Uniprot            | P30443  |
| GeneID             | 3105  |
| Calculated MW      | 38  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only