

ADRB3 Polyclonal Antibody PE Conjugated

Catalog No: #C90557PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADRB3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	310-408/408
Conjugates	PE
Target Name	ADRB3
Other Names	ADRB3
Accession No.	Uniprot:P25962Gene ID:155
Uniprot	P25962
GeneID	155
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only