

Cardiac Troponin T Polyclonal Antibody PE Conjugated



Catalog No: #C90559PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cardiac Troponin T Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-298/298
Conjugates	PE
Target Name	Cardiac Troponin T
Other Names	Cardiac Troponin T
Accession No.	Uniprot:P45379Gene ID:7139
Uniprot	P45379
GeneID	7139
Calculated MW	36
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only