

## Podoplanin Polyclonal Antibody APC Conjugated

Catalog No: #C90566APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Podoplanin Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ms Rat
Immunogen Type	31-130/166
Conjugates	APC
Target Name	Podoplanin
Other Names	Podoplanin
Accession No.	Uniprot:Q86YL7Gene ID:14726
Uniprot	Q86YL7
GeneID	14726
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only