

## CD7 Polyclonal Antibody FITC Conjugated

Catalog No: #C90567FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD7 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms
Immunogen Type	31-100/240
Conjugates	FITC
Target Name	CD7
Other Names	CD7
Accession No.	Uniprot:P09564Gene ID:924
Uniprot	P09564
GeneID	924
Calculated MW	24
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only