

## NMDAR2C Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90568PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NMDAR2C Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/1236
Conjugates	PE-Cy5
Target Name	NMDAR2C
Other Names	NMDAR2C
Accession No.	Uniprot:Q14957Gene ID:2905
Uniprot	Q14957
GeneID	2905
Calculated MW	132
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only