

CD8 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90571PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD8 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms
Immunogen Type	101-200/247
Conjugates	PE-Cy5
Target Name	CD8
Other Names	CD8
Accession No.	Uniprot:P01731Gene ID:12525
Uniprot	P01731
GeneID	12525
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only