

phospho-EphB1/EphB2 (Tyr594) Polyclonal Antibody FITC Conjugated

Catalog No: #C90585FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-EphB1/EphB2 (Tyr594) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	KI(p-Y)ID
Conjugates	FITC
Target Name	phospho-EphB1/EphB2 (Tyr594)
Other Names	phospho-EphB1/EphB2 (Tyr594)
Accession No.	Uniprot:P54762Gene ID:2047
Uniprot	P54762
GeneID	2047
Calculated MW	110
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only