

Laminin alpha 5 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90592Bio

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Laminin alpha 5 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/3717
Conjugates	Bio
Target Name	Laminin alpha 5
Other Names	Laminin alpha 5
Accession No.	Uniprot:O15230Gene ID:3911
Uniprot	O15230
GeneID	3911
Calculated MW	402
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only