

Cystatin-C Polyclonal Antibody PE Conjugated

Catalog No: #C90593PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin-C Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Conjugates	PE
Target Name	Cystatin-C
Other Names	Cystatin-C
Accession No.	Uniprot:P01034Gene ID:1471
Uniprot	P01034
GeneID	1471
Calculated MW	13
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only