

## LRP1 Polyclonal Antibody PE Conjugated

Catalog No: #C90595PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LRP1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2501-2600/4544
Conjugates	PE
Target Name	LRP1
Other Names	LRP1
Accession No.	Uniprot:Q07954Gene ID:4035
Uniprot	Q07954
GeneID	4035
Calculated MW	85/515
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only