

## Aquaporin 5 Polyclonal Antibody FITC Conjugated

Catalog No: #C90597FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aquaporin 5 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/265
Conjugates	FITC
Target Name	Aquaporin 5
Other Names	Aquaporin 5
Accession No.	Uniprot:P55064Gene ID:362
Uniprot	P55064
GeneID	362
Calculated MW	29
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only