

Clostridium difficile Toxin B Polyclonal Antibody FITC Conjugated



Catalog No: #C90601FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Clostridium difficile Toxin B Polyclonal Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Clostridium difficile |
| Immunogen Type | 301-400/2367 |
| Conjugates | FITC |
| Target Name | Clostridium difficile Toxin B |
| Other Names | Clostridium difficile Toxin B |
| Accession No. | Uniprot:P18177Gene ID:4914074 |
| Uniprot | P18177 |
| GeneID | 4914074 |
| Calculated MW | 260 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only