

Actin Polyclonal Antibody Cy5.5 Conjugated

Catalog No: #C90602Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Actin Polyclonal Antibody Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms
Immunogen Type	3-100/377
Conjugates	Cy5.5
Target Name	Actin
Other Names	Actin
Accession No.	Uniprot:P68133Gene ID:58
Uniprot	P68133
GeneID	58
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only