

## HRH4 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90606Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HRH4 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Ms Rat
Immunogen Type	101-200/391
Conjugates	Bio
Target Name	HRH4
Other Names	HRH4
Accession No.	Uniprot:Q91ZY2Gene ID:225192
Uniprot	Q91ZY2
GeneID	225192
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only