

PD-L1 Polyclonal Antibody FITC Conjugated

Catalog No: #C90612FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PD-L1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Ms Rat
Immunogen Type	201-290/290
Conjugates	FITC
Target Name	PD-L1
Other Names	PD-L1
Accession No.	Uniprot:Q9EP73Gene ID:60533
Uniprot	Q9EP73
GeneID	60533
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only