

## Cytokeratin 8 Polyclonal Antibody APC Conjugated

Catalog No: #C90617APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytokeratin 8 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/483
Conjugates	APC
Target Name	Cytokeratin 8
Other Names	Cytokeratin 8
Accession No.	Uniprot:P05787Gene ID:3856
Uniprot	P05787
GeneID	3856
Calculated MW	53
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only