

PTHLH Polyclonal Antibody HRP Conjugated

Catalog No: #C90618HRP



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PTHLH Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	101-177/177
Conjugates	HRP
Target Name	PTHLH
Other Names	PTHLH
Accession No.	Uniprot:P12272Gene ID:5744
Uniprot	P12272
GeneID	5744
Calculated MW	16
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only