

PCDH12 Polyclonal Antibody FITC Conjugated

Catalog No: #C90623FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCDH12 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/1184
Conjugates	FITC
Target Name	PCDH12
Other Names	PCDH12
Accession No.	Uniprot:Q9NPG4Gene ID:51294
Uniprot	Q9NPG4
GeneID	51294
Calculated MW	126
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only