

PCDH18 Polyclonal Antibody PE Conjugated

Catalog No: #C90625PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCDH18 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	381-480/1135
Conjugates	PE
Target Name	PCDH18
Other Names	PCDH18
Accession No.	Uniprot:Q9HCL0Gene ID:54510
Uniprot	Q9HCL0
GeneID	54510
Calculated MW	123
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only