

Oct-3/Oct-4 Polyclonal Antibody PE Conjugated

Catalog No: #C90629PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oct-3/Oct-4 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/360
Conjugates	PE
Target Name	Oct-3/Oct-4
Other Names	Oct-3/Oct-4
Accession No.	Uniprot:P20263Gene ID:5460
Uniprot	P20263
GeneID	5460
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only