

glypican 3 Polyclonal Antibody Bio Conjugated

Catalog No: #C90633Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	glypican 3 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	glypican 3
Other Names	glypican 3
Accession No.	Uniprot:P51654Gene ID:2719
Uniprot	P51654
GeneID	2719
Calculated MW	65
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only