

## AMPK alpha 1 Polyclonal Antibody FITC Conjugated

Catalog No: #C90636FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AMPK alpha 1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein G
Applications	FC
Species Reactivity	Hu Ms Rat
Immunogen Type	351-450/559
Conjugates	FITC
Target Name	AMPK alpha 1
Other Names	AMPK alpha 1
Accession No.	Uniprot:Q13131Gene ID:5562
Uniprot	Q13131
GeneID	5562
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only