

VGLUT1/BNP1 Polyclonal Antibody FITC Conjugated

Catalog No: #C90637FITC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	VGLUT1/BNP1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/560
Conjugates	FITC
Target Name	VGLUT1/BNP1
Other Names	VGLUT1/BNP1
Accession No.	Uniprot:Q9P2U7Gene ID:57030
Uniprot	Q9P2U7
GeneID	57030
Calculated MW	62
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only