

## GPR49 Polyclonal Antibody FITC Conjugated

Catalog No: #C90640FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR49 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	810-907/907
Conjugates	FITC
Target Name	GPR49
Other Names	GPR49
Accession No.	Uniprot:Q9Z1P4Gene ID:8549
Uniprot	Q9Z1P4
GeneID	8549
Calculated MW	98
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only