

EMR1/ F4/80 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90642PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EMR1/ F4/80 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/886
Conjugates	PE-Cy5.5
Target Name	EMR1/ F4/80
Other Names	EMR1/ F4/80
Accession No.	Uniprot:Q14246Gene ID:2015
Uniprot	Q14246
GeneID	2015
Calculated MW	95
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only