

## MRF/C11orf9 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90645Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MRF/C11orf9 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/1151
Conjugates	Bio
Target Name	MRF/C11orf9
Other Names	MRF/C11orf9
Accession No.	Uniprot:Q9Y2G1Gene ID:745
Uniprot	Q9Y2G1
GeneID	745
Calculated MW	124
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only